

PADRE CONCEICAO COLLEGE OF ENGINEERING


EXPERIMENT NO: 14

DETERMINATION OF HARDNESS OF WATER

AIM: To determine the hardness of a given sample of water using EDTA method.

REQUIREMENTS: Eriochrome black-T indicator, 0.01N ZnSO_4 solution, EDTA solution, buffer and water sample

EQUATION:



PROCEDURE: Fill the burette with 0.01N EDTA solution. Pipette 25ml of ZnSO_4 in the conical flask, add 3ml of buffer solution and 3 drops of indicator, then titrate it against EDTA solution till the colour changes from wine red to light blue. Note down the burette reading. Take three such readings.

OBSERVATION:

Burette solution : 0.01N EDTA solution
Pipette solution : 0.01N ZnSO_4 solution
Indicator : Eriochrome black-T
Colour change : wine red to light blue

1 st Reading in 'ml'	2 nd Reading in 'ml'	3 rd Reading in 'ml'	Mean in 'ml'
24.9	24.8	24.9	

CALCULATIONS:

$$N_1 V_1 = N_2 V_2$$

Where N_1 = Normality of EDTA solution

V_1 = Volume of EDTA solution

N_2 = Normality of $ZnSO_4$ solution

V_2 = Volume of $ZnSO_4$ solution

The normality of EDTA solution = _____ N

Part II : Pipette 25ml of water sample in a conical flask, add 3 ml of buffer solution, 2-3 drops of indicator and titrate this against EDTA solution from the burette till the colour changes from wine red to blue, repeat the procedure thrice. Let the mean reading be 'X' ml.

OBSERVATIONS:

Burette solution : 0.01N EDTA solution

Pipette solution : 25ml water

Indicator : Eriochrome black-T

Colour change : wine red to blue

1 st Reading in 'ml'	2 nd Reading in 'ml'	3 rd Reading in 'ml'	Mean in 'ml'
1.9	1.9	1.8	

CALCULATIONS

1000 ml of 1M EDTA = 100gm of $CaCO_3$

1ml of 0.01M EDTA = 1mg of $CaCO_3$

But 1mg/l of $CaCO_3$ = 1 ppm of $CaCO_3$

1ml of 0.01M EDTA = 1 ppm of $CaCO_3$

'X' ml of N_1 EDTA = $\frac{1 \times X \times N_1}{0.01}$

= _____ ppm.

RESULT: The hardness of the given water sample is _____ ppm.